

Kerry
BioHaze™
Haze for days
in your IPAs



Leading to Better

kerry
.com

BioHaze is created from completely natural materials and delivers a cloudy appearance in beer. It acts by binding with protein molecules in beer and forming polyphenol-protein complexes, providing a stable cloudy appearance.

Application & Usage Rate

Recommended dose rate: 60-120 g/bbl (50-100 g/hl, or 500 to 1000 ppm). Equal to 2.5 l of a 2% solution.

Preparation of a 2% Solution:

BioHaze disperses easily in cold water – stir thoroughly to ensure complete dissolution

- Place 264 gal (1000 l) of cold water in tank and turn on the agitator until a good vortex is formed.
- Slowly add 20 kg of BioHaze at the vortex
- Stir until dissolved



Select Ingredients

800.374.2739
bsgcraft.com
sales@bsgcraft.com

© Kerry 2020